

Disposable Sterile Surgical Drape

一、Product name, model and specification

Product name: Disposable sterile surgical drape

Model: The use of the product in the operating room can be divided into: type A, B, C, D, E, F, H, I; The raw materials and process processing classification is divided into: type I, type II, type III and type IV.

Specification: Length(200mm~3000mm)×Width(200mm~3000mm)(See the product minimum packaging label for the specific specifications)

二、Name, address, production address, contact information and after-sales service unit of the registered person/manufacturer

1、Name: Hunan Zhenheyikang Medical Instrument Co.,LTD

2、Address: No.6 Building Jingxiang Energy, No.55 Xianguang east road, Gaoxin district Xiangtan, Hunan

3、Production address: NO.6 Building, NO.8 building, 3rd and 4th floor Jingxiang Energy, No.55 Xianguang east road, Gaoxin district Xiangtan, Hunan, NO 19 Chutian Road, Yuetang District, Xiangtan City, Hunan

4、Contact number: +86 17711627084

5、After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、Manufacturing License NO: Hunan Medical Products Administration Production Licence No. 20150013

四、Product Registration Certificate No, Product Technical Requirements No: HNFDA Certified 20192140275

五、Structure and composition, scope of application

The product is made up of nonwovens or water absorbent layer and liquid resistance layer. The water absorption layer is toilet paper or wood pulp cotton material, and the liquid resistance layer is PE film material. The nonwovens meet the relevant requirements of FZ / T64004 thin adhesive legal nonwovens, and the toilet paper meets the relevant requirements of GB20810-2006 toilet paper (including toilet paper base paper). The PE film meets the relevant requirements of GB13735-2017 polyethylene blown molding agricultural ground covering film.

For clinical cross-contamination.

六、Use method: open the bag can be used

七、Contraindications, precautions, warnings and reminders

Product contraindications: None

Damaged packaging prohibited use, expired products prohibited use;

Disposable use, please destroy after use

Avoid high temperature, humidity, direct sunlight

八、Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用



环氧乙烷灭菌



无菌

九、Installation and use instructions or illustrations

This product does not need to be installed, open the bag can be used.

十、Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in relative humidity $\leq 80\%$, non-corrosive gas and well-ventilated indoor, ordinary transport, in the transportation, storage process should be kept clean and pollution-free, strictly prohibited rain and snow immersion.

十一、Date of production, shelf life or expiration date

Valid for two years, production batch number, production date, expiry date please refer to the small package or qualification label.

十二、 Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description

This product has n o accessories

十三、 Date and version of preparation or revision of the specification

This specification was revised in February 2023, version 2023.