

Medical Nursing Pad

一、Product name, model and specification

Product name : Medical nursing pad

Model : Flat shop type、Paste type、active type、S type、U type、triangular form

Specification : 160*400、210*490、220*500、235*580、240*570、250*330、270*360、270*720、280*400、280*540、280*600、300*400、300*450、300*700、300*780、320*440、320*450、320*460、320*490、320*520、320*540、320*560、320*620、320*630、320*820、330*250、330*450、360*650、380*690、500*400、560*560、600*400、600*450、600*600、600*700、600*750、600*800、600*900、600*1200、600*1500、650*750、650*800、700*1100、700*1200、700*1500、700*1800、750*650、750*900、750*920、750*1000、750*1500、750*2000、800*650、800*750、800*900、800*960、800*1000、800*1200、800*1500、800*1600、820*1015、840*1030、880*750、900*900、900*1200、900*1950、1000*1500、1000*2000、1000*2300、1000*2500.heterotype pad length 500-900 (Unit : mm)。

二、Name、address、production address、contact infomation and after-sales service unit of the registered person/manufacturer

1、Name: Hunan Zhenheyikang Medical Instrument Co.,LTD

2、Address: No.6 Building Jingxiang Energy, No.55 Xianguang east road, Gaoxin district Xiangtan, Hunan

3、Production address:NO.6 Building, NO.8 Building,5th floor Jingxiang Energy, No.55 Xianguang east road,Gaoxin district Xiangtan, Hunan

4、Contact number: +86 17711627084

5、After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、Production Registration Certificate No: Xiangtan Medical Products Administration Production Registration NO 20140004

四、Product Registration Certificate No、Product Technical Requirements No: XTFDA Certified 20230017

五、Product description, intended use

Made of nonwoven fabric and plastic film. Non-sterile provided, single use.

Sanitary care products used in the hospital bed or examination bed.

六、Use method: Open the bag can be used, directly to treat the wound

七、Contraindications, precautions, warnings and reminders

Product contraindications: Please stop using if allergic to the product material, skin surface damaged or part skin verification infected.

Avoid high temperature, humidity, direct sunlight

Damaged packaging prohibited use, expired products prohibited use

Disposable use, please destroy after use

八、Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用

九、Installation and use instructions or illustrations

Choose the appropriate model and specification and use them directly.

十、Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in a non-corrosive gas and well-ventilated room, transported by ordinary transportation means, keep clean and pollution-free, rain and snow immersion.

十一、Date of production, shelf life or expiration date

Valid for three years, production batch number, production date, expiry date please refer to the small package or qualification label.

十二、Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description

This product has n o accessories

十三、Date and version of preparation or revision of the specification

This specification was revised in March 2023, version 2023.