



中国认可
国际互认
检测
TESTING
CNAS L0608

天纺标检测认证股份有限公司

TianFangBiao Standardization Certification & Testing Co.,Ltd.

国家服装质量监督检验中心 (天津)

National Clothing Quality Inspection & Supervision Center(Tianjin)

国家针织产品质量监督检验中心

National Knitted Product Quality Inspection & Supervision Center


检验检测报告

Test Report



T T T S - W T 2 0 2 4 5 0 8 9

Information Provided by Client	Applicant	HUNAN ZHENHEYIKANG MEDICAL INSTRUMENT CO.,LTD		Contact	/	
		NO.6 Building Jingxiang Energy,NO.55 Xiaguang East Road,Gaoxin District, Xiangtan, Hunan		Tel.	/	
	Manufacturer	HUNAN ZHENHEYIKANG MEDICAL INSTRUMENT CO.,LTD NO.6 Building Jingxiang Energy,NO.55 Xiaguang East Road,Gaoxin District, Xiangtan, Hunan				
	Information of Submitted Sample	Sample Name	Medical Surgical Mask	Trademark	/	
		Sample Count	50Piece(s)			
Size		Ear-loop Type	Colour	/		
	Quality Grade	/	Safety Category	/		
	Style No. or Order No.	/				
Test Standards	ASTM F2100-19 Standard Specification for Performance of Materials Used in Medical Face Masks(Level 3 Barrier)					
Test Part Description	1# Blue Mask					
Test Type	Commission Test	Date of Submission	2020-09-18	Date of Completion	2020-09-22	
Test Date	2020-09-18		To	2020-09-21		
Test Standards	See next page(s)					
Conclusion	Test results and compliance refer to next page(s). Stamp of Inspection Unit					
Remarks	Signature Field of He Zhen:Microbe test item;Signature Field of Shan Xuelei:All test items except for Microbe test item.					



Approver

李学雷

Checker

吴春英

Editor

王



中国认可
国际互认
检测
TESTING
CNAS L0608

天纺标检测认证股份有限公司

TianFangBiao Standardization Certification & Testing Co.,Ltd.

国家服装质量监督检验中心 (天津)

National Clothing Quality Inspection & Supervision Center(Tianjin)

国家针织产品质量质量监督检验中心

National Knitted Product Quality Inspection & Supervision Center

检验检测报告

Test Report

TTTS-WT20245089

Page 2 of 4

Test Items	Description	Unit	Standard Requirement	Results	Conclusions	Test Method/Remarks
1# Blue Mask						
Differential Pressure	/	mm H ₂ O/cm ²	<6.0	4.1	Pass	EN 14683:2019+AC:2019 (E)
Resistance to Penetration by Synthetic Blood	/	mmHg	≥160	≥160	Pass	ASTM F1862/F1862M-17
Bacterial Filtration Efficiency	/	%	≥98	1: 99.92 2: 99.79 3: 99.88 4: 99.92 5: 99.79	Pass	ASTM F2101-19
Flammability	/	/	Class 1	Class 1	Pass	US CPSC 16 CFR PART 1610-2019



中国认可
国际互认
检测
TESTING
CNAS L0608

天纺标检测认证股份有限公司
TianFangBiao Standardization Certification & Testing Co.,Ltd.

国家服装质量监督检验中心 (天津)
National Clothing Quality Inspection & Supervision Center(Tianjin)
国家针织产品质量质量监督检验中心
National Knitted Product Quality Inspection & Supervision Center

检验检测报告

Test Report

TTTS-WT20245089

Page 3 of 4

Sample





中国认可
国际互认
检测
TESTING
CNAS L0608

天纺标检测认证股份有限公司

TianFangBiao Standardization Certification & Testing Co.,Ltd.

国家服装质量监督检验中心 (天津)

National Clothing Quality Inspection & Supervision Center(Tianjin)

国家针织产品质量监督检验中心

National Knitted Product Quality Inspection & Supervision Center

检验检测报告

Test Report

【注意事项】

POINTS FOR ATTENTION

Page 4 of 4


- 1.检验报告无“检验专用章”无效。
Report is invalid without stamp of "special seal for the test report".
- 2.复制报告未重新加盖“检验专用章”无效。
Copy report is invalid without re-stamp of "special seal for the test report".
- 3.报告无编写、审核、批准人签字无效。
Report is invalid without collective signatures by editor, checker and approver.
- 4.检验报告涂改无效。
Report is invalid if altered.
- 5.检验报告或复制报告未加盖骑缝章无效 (报告页数多于1页时)。
Report and copy report is invalid without stamp of "Paging Seal" (When the page number more than 1).
- 6.委托检验仅对来样负责,不承担其他连带责任。
Unless otherwise stated the results shown in this report refer only the sample(s) tested.
- 7.对于检验结果若有异议,应于收到报告之日起十五日内向本机构提出,逾期不予受理。
Objection should be issued in 15 days upon receiving the report, overdue opinion is inadmissible.
- 8.未经本机构书面批准,部分复制报告无效。
Part copy report is invalid without the approval of the written documents of the testing organization.

-----报告结束 End of report-----

注意事项以中文为准 The English edition is for reference only

天纺标集团检测单位与地址 Tianfangbiao Groups Others Testing Location

天纺标检测认证股份有限公司
Tianfangbiao Standardization Certification & Testing Co., Ltd.
国家服装质量监督检验中心 (天津)
China National Clothing Quality Inspection & Supervision Center (Tianjin)
国家针织产品质量监督检验中心
China National Knitted Product Quality Supervision Testing Center
地址:天津市南开区鹊桥路25号
Address: No.25, Queqiao Road, Nankai District, Tianjin, China



扫码下载电子报告



扫码关注天纺标