

Disposable Medical mask

一、Product name, model and specification

Product name: Disposable medical mask

Model specification: Ear hanging type: 75X125、90X145、95X145、100X180, Branding type: 100X180 (Unit mm)

二、Name of the registrant / production enterprise, domicile / production address, contact information and after-sales service unit

1、Name: Hunan Zhenheyikang Medical Instrument Co.,LTD

2、Address: No.6 Building Jingxiang Energy, No.55 Xiaguang east road, Gaoxin district Xiangtan, Hunan

3、Production address: NO.6 Building, NO.8 building, 3rd and 4th floor Jingxiang Energy, No.55 Xiaguang east road, Gaoxin district Xiangtan, Hunan, NO 19 Chutian Road, Yuetang District, Xiangtan City, Hunan

4、Contact number: +86 17711627084

5、After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、Manufacturing License NO: Hunan Medical Products Administration Production Licence No. 20150013

四、Product Registration Certificate No, Product Technical Requirements No: HNFDA Certified 20182140010

五、Structure and composition, scope of application

The product is mainly composed of mask body, nose clip and mask belt. The mask body is divided into three layers: inner, middle and outer. The inner layer and outer layer are polypropylene spin-bonded non-woven fabric, the middle layer is polypropylene melt-blown non-woven fabric, and the nose clip is composed of deformable material.

Ethylene oxide sterilization, sterile supply.

Used to wear, block the oral and nasal exhalation or discharge of pollutants in the general medical environment.

六、Use method: open the bag can be used

七、Contraindications, precautions, warnings and reminders

Product contraindications: None

Damaged packaging prohibited use, expired products prohibited use;

Disposable use, please destroy after use

Avoid high temperature, humidity, direct sunlight

八、Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用



环氧乙烷灭菌



无菌

九、Installation and use instructions or illustrations

Disposable medical masks have three layers, and there are inside and outside. The light color surface has moisture absorption function, and should be attached to the mouth and nose, so that the dark face faces outward. With the nose clip facing up, remember not to wear it backwards. According to your face shape, expand the folding surface fully, completely wrap the nose, mouth and jaw, and then press the nose clip to fully fit the mask to the face.

十、Product maintenance and maintenance methods, special storage, transportation conditions, methods

Products should be stored in relative humidity $\leq 80\%$, non-corrosive gas and well-ventilated indoor, ordinary transport, in the transportation, storage process should be kept clean and pollution-free, strictly prohibited rain and snow immersion.

十一、 Date of production, shelf life or expiration date

Valid for two years, production batch number, production date, expiry date please refer to the small package or qualification label.

十二、 Accessories list, including accessories, appendix, lossitage replacement cycle and replacement method description

This product has no accessories

十三、 Date and version of preparation or revision of the specification

This specification was revised in February 2023, version 2023.